



Donor lounge Comfort

Economic and Reliable. Developed to give comfort and security during long treatments such as Dialysis, Blood transplantation, Oncology and other medical treatments.

Technical data

Power supply: 110-230 V 50-60 Hz
 Protection category: IP44
 Handset and MCL protection category: IP66

Standard Dimensions

Seat width: 54 cm
 Total width including arm rests: 88 cm
 Height adjustable from 55 to 75 cm
 Total standard length: 205 cm
 Total standard optional length: 217 cm

Logistic Data Net

Weight: 80 kg
 Gross weight with pallet and box: 120 kg

Features

- Easy to use and maintain
- New design epoxy metal frame equipped with 3 motors adjustable head cushion Handset
- Central brake and emergency button to increase patient's security
- Trendelenburg position - 12 ° operated on the handset
- Manual foot rest adjustment
- Back handle to move the chair easily and hook for the cable
- Electric height adjustment from 55 cm to 75 cm
- Electric back rest adjustment
- Electric leg rest adjustment
- Steel frame structure with light grey epoxy finish
- Multi position arm rests
- Arm rest inclination adjustment
- Upholstery with seams, 12 cm foam thickness
- 4 x 10 cm castors with individual locking system
- Tested and approved patient weight capacity: 250 kg (tested at a dynamic weight of 250 kg and a static weight of 500 kg)

Advantages

Comfort

Foam thickness: 12 cm
 3 electrical motors: back rest, variable height and leg rest
 Smooth electrical motions

Accessibility

Variable height from 55 to 75 cm

Security

Red key on the handset to operate the Trendelenburg position
 Granted stability

Hygiene

Coated fabric Upholstery
 Easy to clean and maintain

Available colors



Options

- Arm rest support: stainless steel
- Foot rest: with metal structure protection and plastic protection options
- Extra comfort upholstery
- Coated fabric M2 or M1
- Back movement handle
- Paper roll holder
- IV serum pole
- Battery
- Crash button
- Foot rest motor
- Foldable ABS tray