




Donor lounge Comfort II

Light and flexible donor lounge. Developed to give comfort and security during long treatments such as Dialysis, Blood transplantation, Oncology and other medical treatments.


Available colors



OLIVE 01173047 CANARD 01173046 MYOSOTIS 01173033 SAPHIR 011 73 005



SABLE 01173008 PEPPER 01173044 BROWN 01173043 NEIGE 01173015



MACADAM 01173007 BRUME 01173006 LIGHT GREY 01173045 POTERIE 01173034



CERISE 01173009 GRENAT 01173002 CHARBON 01173001

Features

- Easy to use and maintain
- New design epoxy metal frame equipped with 3 motors, adjustable head cushion, handset
- Trendelenburg position - 12 ° operated on the handset
- Manual foot rest adjustment
- Electric height adjustment from 55 cm to 75 cm
- Electric back rest adjustment
- Electric leg rest adjustment
- Steel frame structure with light grey epoxy finish
- Multi position arm rests
- Arm rest inclination adjustment
- Upholstery with seams, 12 cm foam thickness
- 4 x 10 cm castors with individual locking system
- Tested and approved patient weight capacity: 250 kg (tested at a dynamic weight of 250 kg and a static weight of 500 kg)

Advantage

Comfort

Foam thickness : 12 cm
3 electrical motors: back rest, variable height and leg rest
Smooth electrical motions

Accessibility

Variable height
from 55 to 75 cm

Hygiene

Coated fabric Upholstery
Easy to clean and maintain

Pricing

Good price/quality relationship

Options

IV serum pole
Foot rest motor

Security

Brake on each individual wheel
Granted stability

Technical Specifications

Power supply: 110-230 V 50-60 Hz
Protection category: Ip44
Handset and MCL protection category: IP66

Logistic Data Net

Weight: 75 Kg
Gross weight: 110 Kg

Standard Dimensions

Seat width: 55 cm
Overall width: 88 cm
Height: 55 - 75 cm
Length: min 198 - max 215 cm

Lmb Technologie GmbH
Möslstr. 17
D-85445 Schwaig

Tel: +49-8122-880 96-0
Fax: +49-8122-880 96-60

www.lmb.de
info@lmb.de