

# Mobile donor lounge DL100



The perfect lounge ...

... for every blood donation center

## Features

- **Portable:** The donor lounge is foldable.
- **Shock position:** A collapsed donor can be set into a shock position supported by a gas tube.
- **Donor slip:** Easy to clean and disinfect. The material has a net structure and is comfortable for the donor.
- **Head rest:** The upholstery of the lounge is integrated in the unit.
- **Arm rest:** The armrest is adjustable on both sides.



Easy to carry on rolls



Adjustable armrest

**Shock position**  
Gas pressure supported



## Technical data

<b>Material:</b>	Stainless steel (round pipe) to prevent scratches and accidents
<b>Weight without arm rest:</b>	17,8 kg
<b>Arm rest:</b>	1,7 kg
<b>Max. net load:</b>	130 kg
<b>Dimensions:</b>	125 x 55 x 22 cm (HxWxD)
<b>High of seat:</b>	58 cm
<b>High of arm rest:</b>	75 cm
<b>Seat extent:</b>	55 cm
<b>Angle of seating:</b>	approximately 140°