

# Blue Box



## Advantages

- No expensive cooling logistic necessary.
- Keeping of optimal temperature.
- Reliability for electronic power outage and defects.
- Reusable.
- Stackable.
- Sealable.







## 3 L cooling elements available as:

- 4°C accu, fluid blue, fluid-crystal-thermometer 0°C-12°C, shows the actual temperature of the element. For the transport of: red cells (Art.N°: KÜ-10045)
- 22°C accu, fluid green, fluid-cristal-termometer 10°C-40°C, shows the actual temperature of the element For the transport of: platelets, whole-blood, cord-blood (Art.N°: KÜ-10023)
- - 30°C accu, fluid yellow, for the transport of: fresh frozen plasma

## Position of 3 L cooling element:

4 positions above the bloodbags are possible



## Technical

<b>Capacity</b>	30 liter
<b>Outer volume (mm)</b>	464 x 368 x 340
<b>Inner volume (mm)</b>	470 x 285 x 230
<b>Material</b>	Expanded

**Distributed by**  
**Lmb Technologie GmbH**  
 Möslstr. 17  
 D-85445 Schwaig  
 Tel: +49-8122-88096-0  
 Fax: +49-8122-88096-60  
 www.lmb.de  
 info@lmb.de

Produced by