

# Blood Transportation Bag – BTB



Designed for safe inhouse transport of Blood and Blood Components

## Features

- The BTB has a perforated design and can be opened by authorized health staff at the final destination with no need for special equipment. Once opened, the BTB cannot be re-sealed!
- The BTB is provided with double-side tape to ensure the safety of blood transportation from the point of origin to its last destination.
- The BTB is produced in 3 layers: a layer of polystyrene foam wrapped with 2 layers of low density polyethylene to ensure mechanical safety and temperature control.
- The complementary information on the front side of the BTB, when correctly filled, will increase the safety of blood transfusion.
- Safety instructions are provided on the rear side of the BTB.
- For longer periods of transportation, please make sure that the BTB is placed inside an insulated and temperature-controlled box.
- Maximum capacity: 1000 ml

## Optional (minimum order quantity 10.000 pieces)

- Upon country and customer needs the information on the front side of the BTB can be modified and company logo can be printed on the BTB.

Distributed by  
**Lmb Technologie GmbH**

Möslstr. 17  
D-85445 Schwaig  
Tel: +49 - 8122 - 880 96 - 0  
Fax: +49 - 8122 - 880 96 - 60  
www.lmb.de  
info@lmb.de

Produced by

