

Donor lounge Comfort III

Economic and Reliable. Developed to give comfort and security during long treatments such as Dialysis, Blood transplantation, Oncology and other medical treatments.



Features

- Donor Lounge Comfort III uses 40% less space than standard dialysis chair beds
- Easy to use and maintain
- Three main parts: back rest, foot rest and seat
- Back, foot, trendelenburg, up and down movements done by 4 DC motors with hand control
- Depending on patient position, arms can open left and right sides. When you move the chair back, adjustable arms does not disturb the parallel position to the ground
- Foot rest adjustable to patient height
- Brakes on each individual wheel, for easy moving and fixing
- Movable head part, adjustable for patient position
- Steel construction has galvanized and electrostatic painting. All parts are coated with cleanable antibacterial ABS plastic

Technical data

Power supply:	100-240 VAC, 50/60 Hz
Protection class:	I
Protection category:	IPX4
In compliance with:	EN 60601-1

Standard dimensions

Chair Width	
Including Armrests:	88 cm
Length:	177 cm
Back segment angle:	0-71°
Back segment length:	63 cm
Height adjustment	55-75 cm
Seat length	48 cm
Head rest length:	22.5 cm
Head rest adjustment:	15 cm
Leg segment length:	45 cm
Foot rest adjustment:	15 cm
Leg segment angle:	0-35°
Safe maximum load:	170 kg

Logistic data

Weight:	80 kg
Gross Weight:	105 kg
Package dimensions:	161 x 91 x 90 cm

Optional pillows



Available colours



Lmb Technologie GmbH
Möslstr. 17
D-85445 Schwaig

Tel: +49-8122-880 96-0
Fax: +49-8122-880 96-60

www.lmb.de
info@lmb.de

Accessories

- **Standard:** IV Pole
- **Optional:** Battery, Eating table, Head pillow, Pillows