

# Tube stripper

## Automatic Tube stripper TS08



### Features

- New and redesigned tube holder for better tube guidance
- Designed for right and left hand use
- Low weight for easy handling
- Swift and complete stripping of blood bag tubes
- Ideal for all commercial blood bag systems on the market
- Automatically centered tubing fed through pair of rollers by a guidance system
- Preventing coagulation of blood in tubing prior to preparing segments
- Power supply with stabilized output, low ripple & low interference, Short circuit & overload protection
- High efficiency and low energy consumption

### General technical data

<b>Overall dimensions</b> (handle closed)	75 x 57 x 195 mm
<b>Overall dimensions</b> (handle open)	110 x 57 x 195 mm
<b>Max stripping tube diameter</b>	6.5 mm
<b>Stripping wheel diameter</b>	33 mm
<b>Stripping wheels angle opening</b>	25 deg
<b>Tube stripper body diameter</b>	40 mm
<b>Max stripping speed</b>	400 mm/sec
<b>Weight without power supply</b>	~800 g
<b>Operating temperature</b>	+10 – +40 °C
<b>Storage and transport temperature</b>	-20 – +50 °C
<b>Power supply voltage</b>	100-240 V AC
<b>Input frequency</b>	50 Hz / 60 Hz
<b>DC motor supply voltage</b>	12V

